SEPARATION REPORT

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- 1 Analysis of middle molecular peptides in normal and uremic body fluids by high-performance gel chromatography, H.Mabuchi and H. Nakahashi, J. Chromatogr., 224(1981)322
- 2 Estimation of molecular weights of peptides by determination of height equivalent to a theoretical plate in size-exclusion chromatography, Y.Shioya, H.Yoshida and T.Nakajima, J.Chromatogr., 240(1982)341
- 3 High-performance gel-permeation chromatography of polypeptides in a volatile solvent : Rapid resolution and molecular weight estimations of proteins and peptides on a column of TSK-G3000-PW, G.D. Swergold and C.S.Rubin, Anal. Biochem., 131(1983) 295
- 4 Molecular weight determination of peptides by high-performance gel permeation chromatography, W. O.Richter, B.Jacob and P.Schwandt, Anal. Biochem., 133(1983)288
- 5 Desalting of peptides by high-performance gel permeation chromatography, W.O.Richter and P. Schwandt, J.Chromatogr.,288(1984)212
- 6 Radioimmunoassay for α -human and rat atrial natriuretic polypeptide K.Nakao, A.Sugawara, N. Morii, M.Sadamoto, M.Suda, J.Soneda, T.Ban, M. Kihara, U.Yamori, M.shimokura, Y.Kiso and H. Imura Biochem. Biophys. Res. Comm. 124(1984)815
- 7 General method for the separation of cyanogen bromide digests of proteins by high-performance liquid chromatography ; Rabbit skeletal troponin I, C.T.Mant and R.S.Hodges, J.Chromatogr., 326(1985) 349
- 8 Amino-acid sequence of β -tubulin from Physarum polycephalum myxamoebae, M.S.-Wowra, L. Clayton, P.Dawson, K.Gull and M.Little, Eur. J. Biochem., 161(1986)669
- 9 High performance size exclusion liquid chromatography of small molecular weight peptides from protein hydrolysates using methanol as a mobile phase additive, M.A.Vijayakshmi, L.Lemieux and J. Amiot, J.Liquid Chromatogr.,9(1986)3559
- 10 High-performance gel permeation chromatography of proteins and peptides on columns of TSK-G-2000-SW and TSK-G3000-SW : A volatile solvent giving separation based on charge and size of polypeptides, G.B.Irvine and C.Shaw, Anal. Biochem., 155(1986)141
- 11 Fibrinogens Kawaguchi and Osaka : An amino acid substitution of A α arginine-16 to cysteine which forms an extra interchain disulfide bridge between the two A α chains, T.Miyata, S.Terukina, M.Matsuda, A.Kasamatsu, Y.Takeda, T.Murakami and S. Iwanaga, J.Biochem., 102(1987)93
- 12 Construction of a model for the aggregation and cross-linking region(7S domain) of type IV collagen based upon an evaluation of the primary structure of the $\alpha 1$ and $\alpha 2$ chains in this region, B.Siebold, R. Qian, R.W.Glanville, H.Hofmann, R.Deutzmann and K.Kuhn, Eur. J.Biochem., 168(1987)569
- 13 Size-exclusion high-performance liquid chromatography of peptides : Requirement for peptide standards to monitor column performance and nonideal behaviour, C.T.Mant, J.M.R.Parker and R.S. Hodges, J.Chromatogr., 397(1987)99

- 14 The combined use of size-exclusion and reversed-phase high-performance liquid chromatography for characterization of β -endorphin processing pathways, A.I. Smith, C.A. Wallace, M.C. Cheng and D.J. Autelitano, *J. Chromatogr.*, 416(1987)255
- 15 Functional homology of 6-phosphofructo-2-kinase/fructose - 2,6 - bisphosphatase, phosphoglycerate mutase, and 2,3-bisphosphoglycerate mutase, A. Tauler, M.R. El-Maghrabi and S.J. Pilkis, *J. Biol. Chem.*, 262(1987)16808
- 16 Isolation and analysis of peptidic fragments of α -gliadin using reversed phase high-performance liquid chromatography, P. Kocna, P. Fric, M.K. Holaková, J. Slaby, E. Kasárek and W. Th. J. M. Hekkens, *J. Chromatogr.*, 434(1988)429
- 17 Analytical strategy for determination of active site sequences in aminoacyl-tRNA synthetases, C. Beauvallet, C. Hountoundji and J.-M. Cshmitter, *J. Chromatogr.*, 438(1988)347
- 18 Purification and partial biologic characterization of a human lymphocyte derived peptide with potent neutrophil-stimulating activity, J.-M. Schroder, U. Mrowietz and E. Christophers, *J. Immunol.*, 140(1988)3534
- 19 Amino acids of the Torpedo marmorata acetylcholine receptor α subunit labelled by a photoaffinity ligand for the acetylcholine binding site M. Dennis, J. Giraudat, F.K. Hibert, M. Goeldner, C. Hirth, J.Y. Chang, C. Lazure, M. Chretien and J.-P. Changeux, *Biochemistry*, 27(1988)2346
- 20 Immunogenicity and Allergenicity of whey protein hydrolysates, J. Asselin, J. Amiot, S.F. Gauthier, W. Mourad and J. Hebert, *J. Food Sci.*, 53(1988)1208
- 21 Optimization of peptide separations in high-performance liquid chromatography, C.T. Mant and R.S. Hodges, *J. Liquid Chromatogr.*, 12(1989)139
- 22 Identification of peptides from a peptic haemoglobin hydrolysate produced at pilot-plant scale, by high-performance liquid chromatography and mass spectrometry, J.M. Piot, D. Guillochon, Q. Zhfao, G. Ricart, B. Fournet and D. Thomas, *J. Chromatogr.*, 481(1989)221
- 23 Purification and characterization of human neutrophil peptide 4, a novel member of the defensin family, C.G. Wilde, J.E. Griffith, M.N. Marra, J.L. Snable and R.W. Scott, *J. Biol. Chem.*, 264(1989)11200
- 24 Cross-linking studies of the cholesterol hydroxylation system from bovine adrenocortical mitochondria, I.V. Turko, T.B. Adamovich, N.M. Kirillova, S. A. Usanov and V.L. Chashchin, *Biochim. Biophys. Acta*, 996(1989)37