


# Toyopearl® QAE-550 Resins

(strong anion exchange chromatography resins)


**Part Numbers:**

43271	Toyopearl QAE-550C, 100 $\mu$ m, 100mL	21364	ToyoScreen QAE-650C, 1mL x 6 ea.
14026	Toyopearl QAE-550C, 100 $\mu$ m, 250mL	21365	ToyoScreen QAE-650C, 5mL x 6 ea.
14704	Toyopearl QAE-550C, 100 $\mu$ m, 1L	21392	ToyoScreen IEC Anion Mix Pack, 1mL x 3 grades x 2 ea. (DEAE, SuperQ, QAE)
14027	Toyopearl QAE-550C, 100 $\mu$ m, 5L	21393	ToyoScreen IEC Anion Mix Pack, 5mL x 3 grades x 2 ea. (DEAE, SuperQ, QAE)
18365	Toyopearl QAE-550C, 100 $\mu$ m, 50L	21396	ToyoScreen IEC Mix Pack, 1mL x 6 grades x 1 ea. (DEAE, SuperQ, QAE, CM, SP650, SP-550)
		21397	ToyoScreen IEC Mix Pack, 5mL x 6 grades x 1 ea. (DEAE, SuperQ, QAE, CM, SP650, SP-550)
		21400	ToyoScreen Column Holder
		42194	ToyoScreen Column Holder with fittings

**Product Description**

Toyopearl chromatographic resins are comprised of a methacrylic polymer incorporating high mechanical and chemical stability. Resins are available as non-functionalized "HW" series resins for size exclusion separations, and derivatized with surface chemistries for alternative modes of chromatography such as ion exchange, hydrophobic interaction or affinity separations.

Toyopearl QAE-550 chromatographic resins are designed for anion exchange chromatography. This chromatographic mode separates molecules on the basis of ionic interactions between the sample and the resin. The separation is usually accomplished in buffered aqueous solution with a gradient of increasing ionic strength. Alternatively, pH adjustment may be used for control of elution.

**Formulation**

72% (v/v) slurry in 20% (v/v) ethanol

This sheet contains the suggested operating conditions and specification ranges for Toyopearl QAE-550 resins. Lot-specific specifications are included in the Certificate of Analysis (C of A) shipped with the product.

**A. OPERATING CONDITIONS**

Pressure Drop Across the Column	max. 0.3 MPa (recommended)
Operating Linear Flow Rate	normally 10 – 600 cm/hour
Cleaning in Place/Sanitization	0.1mol/L NaOH or 0.1mol/L HCl
Long Term Storage Conditions	20% (v/v) ethanol (shipping solvent)

**B. SPECIFICATIONS**

Appearance	White resin slurry which settles upon standing
Particle Size Distribution (>80% within range)	50 – 150 $\mu$ m for C-grade
Mean Pore Diameter	50 nm
Exclusion Limit (base resin)	7 x 10 <sup>5</sup> Dalton (globular protein)
Ion Exchange Capacity	0.33 $\pm$ 0.05 eq/L
Protein Adsorption Capacity	70 $\pm$ 10 g/L (Bovine Serum Albumin)
Bacterial Count	$\leq$ 100 CFU/mL
Endotoxin Concentration	$\leq$ 10.0 endotoxin units/mL
Resin Volume per Container	The indicated volume is the settled resin volume.
Foreign Substance	not observed, no visual evidence
Eluable Matter	$\leq$ 0.2%
Shelf Life Stability	Minimum 10 years