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M. Begić et al., Salt-tolerant cation-exchanger containing sulfate groups as a viable alternative for mixed-mode type and heparin-base affinity resin, Biotechnol. J., 2021; 2100100, CC BY NC ND, <https://doi.org/10.1002/biot.202100100>

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FIGURE 1 Difference between sulfonium-type (left) and sulfate-type ion exchanger (right)


FIGURE 2 Heat plot of pH and salt concentration dependence on normalized binding capacity of Lysozyme for: (A) TOYOPEARL GigaCap S-650 M, (B) TOYOPEARL SP-650 M and (C) TOYOPEARL Sulfate-650F

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
FIGURE 4 Purification of clotting factor IX and removal of other vitamin K-dependent clotting factors (FII, FVII, and FX) on sulfate resin TOYOPEARL Sulfate-650F. Conditions: Äkta Explorer chromatographic system (Cytiva) with a 2 mL sample loop, and a 1 mL column, flow 2 mL min⁻¹, room temperature. Buffer A (starting buffer), 400 mM NaCl, 20 mM phosphate, pH 4.5, 39.1 mS cm⁻¹; Collected fractions were eluted with Buffers 1–7 (cf. Materials and methods). Blue: Factor IX activity (IU mL⁻¹) and enrichment (red) were presented. The activity of FX was detected only in the flow-through fraction, and no FX clotting activity was detected in the eluate

Ref.; M. Begić et al., Biotechnol. J., 2021; 2100100, CC BY NC ND, <https://doi.org/10.1002/biot.202100100>

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1. T. Burnouf, An overview of plasma fractionation, Annals of Blood, 2018, 3: 33, <https://aob.amegroups.com/article/view/4496/5240>
2. J. M. Bult, Plasma fractionation technologies: benefit and limitations, September 21, 2021 <https://www.google.com/url?sa=t&rct=j&q=&esrc=s&source=web&cd=&ved=2ahUKEwjV4c3nIbf7AhUQCkHewNDUkQFnoECAsQAQ&url=https://www.isbtweb.org/2Fresource/2Fjan-bult-pdf.html&usq=AOvVaw2kL9acty1TOUXLib1B7sgD>
3. M. Begić et al., Salt-tolerant cation-exchanger containing sulfate groups as a viable alternative for mixed-mode type and heparin-base affinity resin, Biotechnol. J., 2021; 2100100, CC BY NC ND, <https://doi.org/10.1002/biot.202100100>

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