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

Column ; TSKgel G6000PW (7.5 mm I.D. x 30 cm)
 Eluent ; PBS
 Flow rate ; 1.0 mL/min
 Detection ; UV(210 nm), Light scattering
 Temperature ; Room temperature
 Sample ; Erythrocyte Extracellular vesicles

Ref.; Z. Varga et al., Towards traceable size determination of extracellular vesicles, *J. Extracell. Vesicles*, 3 (2014) 23298, CC BY, <https://doi.org/10.3402/jev.v3.23298>

Column ; TSKgel G4000SW_{xL} (7.8 mm I.D. x 30 cm)
 Eluent ; 20 mmol/L phosphate buffer (pH 7.2) containing 0.15 mol/L NaCl
 Flow rate ; 0.5 mL/min
 Detection ; UV(220 nm), Light scattering
 Temperature ; 25 °C
 Sample ; Mesenchymal stem cell (MSC) cultured medium after 24, 48, 72 hr
 フラクションP1はエクソソームと同定された

Ref.; S. S. Tan et al., Therapeutic MSC exosomes are derived from lipid raft microdomains in the plasma membrane, *J. Extracellular Vesicles*, 2 (2013) 22614, CC BY, <https://doi.org/10.3402/jev.v2i0.22614>

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1. Z. Varga et al., Towards traceable size determination of extracellular vesicles, *J. Extracell. Vesicles*, 3 (2014) 23298, CC BY, <https://doi.org/10.3402/jev.v3.23298>
2. G. E. Benedetto et al., Multiple Techniques for Size Determination of Generalized Modules for Membrane Antigens from *Salmonella typhimurium* and *Salmonella enteritidis*, *ACS Omega* 2017, 2, 11, 8282–8289, <https://dx.doi.org/10.1021/acsomega.7b01173>
3. R. M. M. da Silva, A chromatography approach towards large-scale purification of mesenchymal stem/stromal cell-derived extracellular vesicles, (2020) Instituto Superior Técnico, Lisboa
4. S. S. Tan et al., Therapeutic MSC exosomes are derived from lipid raft microdomains in the plasma membrane, *J. Extracellular Vesicles*, 2 (2013) 22614, CC BY, <https://doi.org/10.3402/jev.v2i0.22614>
5. K. Vrabec et al., Emerging tools for exosome: Purification and in-process monitoring, *BioProcess International*, 17(3)E1, March 2019
6. 研究開発の俯瞰報告書、2.4.2 細胞外微粒子・細胞外小胞、ライフサイエンス・臨床医学分野(2021), p499 <https://www.jst.go.jp/crds/pdf/2020/FR/CRDS-FY2020-FR-04.pdf>
7. C. Thery et al., Minimal information for studies of extracellular vesicles 2018 (MISEV2018): a position statement of the International Society for Extracellular Vesicles and update of MISEV2014 guideline, *J. Extracellular Vesicles*, 2018 Vol 7, 1535750, <https://doi.org/10.1080/20013078.2018.1535750>
8. K. Sidhom et al., A review of exosomal isolation methods: is size exclusion chromatography the best option? *Int. J. Mol. Sci.* 2020 Sep 4;21(18):6466, <https://doi.org/10.3390/ijms21186466>
9. H. Kaddour et al., The past, present, and future of the size exclusion chromatography in extracellular vesicles separation, *Viruses*, November 2021, 13(11):2272, <https://dx.doi.org/10.3390/v13112272>

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