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- Center for Biologics Evaluation and Research, Guidance for Industry, Considerations for plasmid DNA vaccines for infectious disease indications, US FDA, November 2007, <https://www.fda.gov/files/vaccines,%20blood%20&%20biologics/published/Guidance-for-Industry--Considerations-for-Plasmid-DNA-Vaccines-for-Infectious-Disease-Indications.pdf>
- J. Whitley et al., Development of mRNA manufacturing for vaccines and therapeutics: mRNA platform requirements and development of a scalable production process to support early phase clinical trials, Volume 242, April 2022, Pages 38-55, <https://doi.org/10.1016/j.trsl.2021.11.009>
- J. Urthaler et al., Improved downstream process for the production of plasmid DNA for gene therapy, Acta Biochimica Polonica, 52 (2005) 703-711, https://doi.org/10.18388/abp.2005_3434
- S. G. L. Quaak et al., Naked Plasmid DNA Formulation: Effect of Different Disaccharides on Stability after lyophilization, AAPS PharmSciTech, 11 (2010) 344-350, <https://doi.org/10.1208/s12249-010-9391-2>
- C. R. Smith et al., Separation of topological forms of plasmid DNA by anion-exchange HPLC; Shifts in elution order of linear DNA, J. Chromatogr. B., 854 (2007) 121-127, <https://doi.org/10.1016/j.jchromb.2007.04.005>
- 東ソー、テクニカルインフォメーション No.205, 陰イオン交換クロマトグラフィーによるプラスミド DNA pBR322 の立体構造異性体の分離 https://www.separations.asia.tosohbioscience.com/File%20Library/TBJS/Lit_JP/TechInfo/TI-205r.pdf
- M. J. Molloy et al., Effective and robust plasmid topology analysis and the subsequent characterization of plasmid isoforms thereby observed, Nucleic Acids Research, 2004 Vol. 32, No. 16, e129, <https://doi.org/10.1093/nar/gnh124>
- 東ソー、セパレーションレポート No. 90, 核酸分離用イオン交換クロマトグラフィー用カラム TSKgel DNA-NPR による DNA フラグメントの分離、https://www.separations.asia.tosohbioscience.com/File%20Library/TBJS/Lit_JP/SepaReport/sr090.pdf
- 東ソー、セパレーションレポート No. 109, 高性能イオン交換クロマトグラフィー用カラム TSKgel STAT シリーズについて、https://www.separations.asia.tosohbioscience.com/File%20Library/TBJS/Lit_JP/SepaReport/SR_No109-STAT-.pdf
- M. Wang et al., Separation of plasmid DNA pBR 322 steric isoforms by HPLC, Nucleic Acids Therapeutics Society of Japan, 2021, poster 100033
- A. Abdulrahman et al., Recent advances in chromatographic purification of plasmid DNA for gene therapy and DNA vaccines: A Review, Analytica Chimica Acta 1025 (2018) 41-57, <https://doi.org/10.1016/j.aca.2018.04.001>

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