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
Ref.; A. Bendelac et al., BioProcess International, October 2022, 20 (10) 23を改変

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1. 東ソー、リーフレット、TSKgel HIC-ADC Butyl.
2. J. Jones et al., Native size-exclusion chromatography-mass spectrometry: suitability for antibody-drug conjugate drug-to-antibody ratio quantitation across a range of chemotypes and drug-loading levels, MABS, 2020, VOL. 12, NO. 1, e1682895, <https://doi.org/10.1080/19420862.2019.1682895>.
3. Y. Matsuda et al., Chromatographic analysis of site-specific antibody-drug conjugates produced by AJICAP first-generation technology using a recombinant FcγIIIa receptor-ligand affinity column, J. Chromatogr. B, 1177 (2021) 122753, <https://doi.org/10.1016/j.jchromb.2021.122753>
4. Y. Fujita et al., Glycoform analysis of intact antibodies using FcR column online mass spectrometer, Tosoh Research & Technology Review, 66 (2022) 47-55, https://www.tosoh.co.jp/technology/assets/22_3_5.pdf
5. A. R. Nanna et al., Generation and validation of structurally defined antibody-siRNA conjugates, Nucleic Acid Research, 48 (10), (2020) 5281-5293, <https://doi.org/10.1093/nar/gkaa286>
6. Y. Matsuda et al., A Purification Strategy Utilizing Hydrophobic Interaction Chromatography to Obtain Homogeneous Species from a Site-Specific Antibody Drug Conjugate Produced by AJICAP™ First Generation, Antibodies 2020, 9 (2), 16, <https://doi.org/10.3390/antib9020016>
7. 東ソー、テクニカルノートトヨパール® No. 6、抗体薬物複合体（ADC）の精製はこれで決まり！
8. 東ソー、テクニカルノートトヨパール® No. 11、抗体薬物複合体（ADC）の精製はこれで決まり！（2）

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