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- 1 Monoclonal antibodies: Purification, fragmentation and application to structural and functional studies of class I MHC antigens, P.Parham, M.J. Androlewicz, F.M.Brodsky, N.J.Holmes and J.P. Ways, *J.Immunol. Methods*, 53(1982)133
- 2 On the fragmentation of monoclonal IgG₁, IgG₂ and IgG_{2b} from BALB/c mice P.Parham, *J.Immunol.*, 131(1983)2895
- 3 Rapid and efficient purification of mouse monoclonal antibodies from ascites fluid using high performance liquid chromatography S.W.Burchiel, J.R. Billman and T.R.Alber, *J.Immunol. Methods*, 69(1984)33
- 4 Single step purification of monoclonal antibody from murine ascites and tissue culture fluids by anion exchange high performance liquid chromatography M.J.Gemski, B.P. Doctor, M.K.Gentry, M.G.Pluskat and M.P.Strickler, *BioTechniques*, 3(1985)378
- 5 One-step procedure for the rapid isolation of mouse monoclonal antibodies and their antigen binding fragments by fast protein liquid chromatography on a Mono Q anion-exchange column, P.Clezardin, J.L. McGregor, M.Manach, H.Boukerche and M.De-chavanne, *J.Chromatogr.*, 319(1985)67
- 6 A High-performance liquid chromatographic procedure for the purification of mouse monoclonal antibodies, J.R.Deschamps, J.E.K.Hildreth, D.Derr and J.T.August, *Anal. Biochem.*, 147(1985)451
- 7 Rapid purification of monoclonal antibodies by high-performance liquid chromatography, B.Pavlu, U. Johansson, C.Nyhlen and A.Wichman, *J. Chromatogr.*, 359(1986)449
- 8 Use of new HPLC resins to solve old problems, D.G. Klapper, S.Osgood, R.Esch and Olson, *J. Liquid Chromatogr.*, 9(1986)1613
- 9 Purification and analysis of monoclonal antibodies from mouse ascites fluid R. Ehrnstrom and B. Gustavsson, *Internat. Lab.*, Sept/Oct(1987)26
- 10 Combination of zeta prep mass ion-exchange media and high-performance cation-exchange chromatography for the purification of high-purity monoclonal antibodies, A. Jungbauer, F.Unterluggauer, F.Steindl, F.Ruker and H.Kattinger, *J.Chromatogr.*, 397(1987)313
- 11 High-performance analytical applications of immobilized metal ion affinity chromatography, M. Belew, T.T. Yip, L.Andersson and R. Ehrnstrom, *Anal. Biochem.*, 164(1987)457
- 12 Purification of ascitic fluid-derived murine monoclonal antibodies by anion-exchange and size-exclusion high-performance liquid chromatography, H. H.Hwang, M.C.Healey and A.V.Johnston, *J. Chromatogr.*, 430(1988)329
- 13 Protein A, hydroxyapatite and diethylaminoethyl: Evaluation of three procedures for the preparative purification of monoclonal antibodies by high-performance liquid chromatography, C.Poiesi, A. Tamanini, S.Ghielmi and A.Albertini, *J. Chromatogr.*, 465(1989)101
- 14 Purification of mouse monoclonal immunoglobulin M by ion-exchange liquid chromatography, S.W.Rothman, M.K.Gentry, R.D.Gawne, A.S.Dobek, R.Ogert, M.J.Stone and M.P.Strickler, *J.Liquid Chromatogr.*, 12(1989)1935
- 15 Separation of monoclonal antibody from ascites fluid by high-performance hydroxyapatite and hydrophobic interaction chromatography, Y. Yamasaki, T.Kitamura and Y.Kato, *The 9th International Symposium on HPLC of Proteins, Peptides and Polynucleotides*, (1989), poster #420

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