

Instructions for TSKgel Boronate-5PW

Support

Boronate-5PW was prepared by introducing m-aminophenylboronic acid into GS000PW.

Particle size : 10 μm

Pore size : ca. 1000 \AA

Column

Size : 75 x 7.5 mm I.D.

Solvent : distilled water

pH range for separation

2 - 9

Salt concentration range

less than 3 M

Flow rate range

less than 1.2 ml/min

Flow rates of 0.5 - 1.0 ml/min are recommended in general.

Column washing

0.1 - 0.2 N NaOH are very effective to wash or regenerate columns.

Usually, columns can be regenerated by injecting 0.1 - 0.2 N NaOH of 1 - 2 ml several times using sample injector. When this procedure did not help, wash the column by injecting 20 - 40 % acetic acid of 1 - 2 ml several times.

Guard column

The use of guard column (10 x 6 mm I.D.) is recommended for a long service life of analytical columns.

The change of both inlet side filter and support once a week is recommended for every day use.

The guard column can be packed by pouring the concentrated slurry of Boronate-5PW guard gel into it and then sucking with aspirator.

The old filter can be re-used after washing by exposing it to ultrasonic in 0.1 - 0.2 N NaOH for about 30 min.

Boronate-5PW guard gel

Support prepared by introducing m-aminophenylboronic acid into GS000PW of 20 - 30 μm in particle diameter for use in guard column

Storage

When Boronate-5PW columns are stored for more than several days, replace the solvent in columns with distilled water.

Store the columns at around 4°C.

Test for theoretical plate number

Following conditions are recommended to test columns for their theoretical plate numbers and As(asymmetry factor).

Eluent : distilled water

Flow rate : 1 ml/min

Sample : acetone (0.5 %, 20 μl)

Boronate-5PW columns have theoretical plate numbers more than 1,300 plates per column and As of 1.0 - 2.0 at a time of delivery.