

3 micron Reversed-phase Chromatography Column:
TSKgel ODS-100V 3 μ m

Table of Contents

	Page
1. Introduction	1
2. Fundamental properties of TSKgel ODS-100V 3 μ m	1
2-1 Properties of the packing materials	1
2-2 H – u curve (Van Deemter curve)	3
2-3 Residual ion exchange activity	5
2-4 Effects of mobile phase on LC/MS(/MS) analysis	7
3. Applications	9
4. Conclusions	12

1. Introduction

Reversed-phase chromatography (RPC) is the most frequently employed separation mode in high-performance liquid chromatography, because it can be used for many different types of sample and has superior operability and separation performance. Reversed-phase octadecyl (C18)-bonded silica gel is widely used as a packing material for examination of low-molecular-weight compounds such as pharmaceuticals.

Tosoh has introduced the TSKgel ODS-100V 3 μ m column, containing packing material with a smaller particle size than that of the previous model, TSKgel ODS-100V 5 μ m (see Separation Report No. 102). To ensure that the new column maintains the separation properties of TSKgel ODS-100V 5 μ m, it contains monolayer octadecyl-bonded (C18) silica gel with a particle size of 3 μ m. TSKgel ODS-100V 3 μ m shares many of the same features as TSKgel ODS-100V 5 μ m, such as superior end-capping efficiency of residual silanol groups, favorable peak shapes for basic and acidic compounds, and the possibility of using mobile phases free from organic solvents. This article discusses the fundamental properties of TSKgel ODS-100V 3 μ m and presents some of applications.

2. Fundamental properties of TSKgel ODS-100V 3 μ m

2-1. Properties of the packing materials

Table 1 compares the fundamental characteristics of TSKgel ODS-100V 3 μ m, TSKgel ODS-100V 5 μ m and TSKgel ODS-100Z 5 μ m. Apart from particle size, the base silica gels of TSKgel ODS-100V 3 μ m and TSKgel

ODS-100V 5 μ m have the same fundamental properties. Also, TSKgel ODS-100V 3 μ m and TSKgel ODS-100V 5 μ m share the same surface modification (functional group bonding and residual silanol group end-capping), and, as shown in **Table 2**, most of the HPLC values are comparable. Furthermore, the carbon content of TSKgel ODS-100V 3 μ m has been adjusted to achieve a retention comparable to that of TSKgel ODS-100V 5 μ m. Therefore, the separation properties of the both columns are comparable.

Figure 1 shows chromatograms of test mixture which contains the same component as NIST SRM870, obtained using TSKgel ODS-100V 3 μ m and TSKgel ODS-100V 5 μ m. The retention of each peak and the asymmetry factors for amitriptyline (a basic compound) and quinizarine (a metal chelating compound) are nearly the same as each other.

Figure 2 shows the relationship between the hydrophobicity of common low-molecular-weight compounds (hydrophobicity parameter: log P) and retention (retention factor: log k') using TSKgel ODS-100V 3 μ m and TSKgel ODS-100Z (3 μ m prototype). Compared to TSKgel ODS-100Z 3 μ m, with high carbon content (20%) and low packing material surface polarity, the retention of compounds with low log P values (hydrophilic compounds) was higher for TSKgel ODS-100V 3 μ m (oval region in the lower left corner). Thus, due to the higher surface polarity of the packing material, hydrophilic compounds are retained more strongly by TSKgel ODS-100V 3 μ m when compared to TSKgel ODS-100Z 3 μ m.

Table 1. Fundamental properties of TSKgel ODS-100V and ODS-100Z

Column	Particle size (μ m)	Pore size (\AA)	Specific surface area (m^2/g)	Pore volume (mL/g)	Functional group	Carbon content ^{*)} (%)	Phase structure
TSKgel ODS-100V 3 μ m	3	100	450	1.10	C18	15	Monolayer
TSKgel ODS-100V 5 μ m	5	100	450	1.10	C18	15	Monolayer
TSKgel ODS-100Z 5 μ m	5	100	450	1.10	C18	20	Monolayer

^{*)} Measured by quantitative elemental analysis
 $\text{\AA} = 1 \times 10^{-10}\text{m}$




Figure 1. Standard sample chromatograms

Column: A) TSKgel ODS-100V 3µm
(4.6 mm I.D. × 15 cm)
B) TSKgel ODS-100V 5µm
(4.6 mm I.D. × 15 cm)

Eluent: 20 mmol/L Phosphate Buffer (pH 7.0)
/CH₃OH = 20 / 80

Flow rate: 1.0 mL/min
Detection: UV 254 nm
Sample: 1. Uracil, 2. Toluene, 3. Ethyl benzene,
4. Quinizarine, 5. Amitriptyline
Inj. volume: 10µL




Figure 2. Relationship between hydrophobicity parameter (log P) and retention (log k')

Column: TSKgel ODS-100V 3µm (4.6 mm I.D. × 15 cm)
TSKgel ODS-100Z (3 µm: prototype)
(4.6 mm I.D. × 15 cm)

Eluent: H₂O/CH₃CN = 60 / 40
Flow rate: 1.0 mL/min
Detection: UV 254 nm
Temperature 40°C
Sample: Aniline, Benzonitorile, Methyl aniline, Anisol,
Methyl benzoate, Benzene, Ethyl benzoate,
Toluene, Ethyl benzene, Xylen

Table 2. HPLC properties of TSKgel ODS-100V and ODS-100Z

Column	Retention coefficient k'	Stereoselectivity α	Hydrogen binding α	Hydrophobicity α	Surface polarity α	Ionization					Coordinate linkage		Retention reduction (RT2Adn/RT1Adn)
						Basic			Acidic		α(k'Quini/k'EB)	As'f(Quini)	
						As'f(Des)	α(k'Ami/k'EB)	As'f(Ami)	α(k'For/k'Ac)	As'f(For)			
TSKgel ODS-100V 3 µm	1.78	1.24	0.47	1.64	0.54	1.62	2.60	1.08	0.48	1.29	1.98	1.02	97.8
TSKgel ODS-100V 5 µm	1.80	1.25	0.45	1.64	0.53	1.59	2.60	1.21	0.48	1.32	1.98	1.16	99.0
TSKgel ODS-100Z 5 µm	2.42	1.31	0.40	1.72	0.43	1.62	2.38	1.07	0.44	1.41	1.77	1.20	-

1. Retention coefficient: k'(Naphthalene)

2. Stereoselectivity: α=k'(Triphenylene)/k'(o-Terphenyl)

3. Hydrogen binding: α= k'(Caffeine)/k'(Phenol)

4. Hydrophobicity: α= k'(Toluene)/k'(Benzene)

5. Surface polarity: α= k'(Methyl benzoate)/k'(toluene)

6. Ionization

1) As'f(Des)=As'f(Desipramine) (pH 7.0)

2) α(k'Ami/k'EB): α=k'(Amitriptyline)/k'(Ethyl Benzene)

3) As'f(Ami)=As'f(Amitriptyline)

4) α(k'For/k'Ac): α=k'(Formic acid)/k'(Acetic acid)

5) As'f(For)=As'f(Formic acid)

7. Coordinate linkage

1) α(k'Quini/k'EB): α=k'(Quinizarine)/k'(Ethyl Benzene)

2) As'f(Quini)=As'f(Quinizarine)

8. Retention: (%)

RT1Adn: Initial elution time for adenine

RT2Adn: Elution time for adenine at 30 minutes after injection

*As: asymmetry coefficient

2-2. H – u curve (Van Deemter curve)

Figure 3 shows the relationship between height equivalent to a theoretical plate (HETP) and linear velocity for TSKgel ODS-100V 3 μ m and TSKgel ODS-100V 5 μ m. Because the particle size of the former is smaller than that of the latter, its HETP is smaller (higher column efficiency). Also, for TSKgel ODS-100V 5 μ m, the smallest HETP and high column efficiency are obtained at a linear velocity of 4 - 6 cm/min, but with TSKgel ODS-100V 3 μ m, which has a smaller particle size, column efficiency is the highest at a greater linear velocity of \geq 6 cm/min. Furthermore, for TSKgel ODS-100V 5 μ m, a greater linear velocity (\geq 6 cm/min) leads to lower column efficiency, while for TSKgel ODS-100V 3 μ m, column efficiency does not decrease much at high linear velocity, and high column efficiency is maintained over a wide range of velocities (6 - 10 cm/min; about 1.0 - 1.7 mL/min for a column with an internal diameter of 4.6 mm). Therefore, compared to TSKgel ODS-100V 5 μ m, chromatography can be performed with higher resolution and higher velocity.




Figure 3. Effects of linear velocity on height equivalent to a theoretical plate (HETP)

Column: TSKgel ODS-100V 3 μ m (4.6 mm I.D. \times 15 cm)
 TSKgel ODS-100V 5 μ m (4.6 mm I.D. \times 15 cm)
 Eluent: H₂O/CH₃OH = 30/70
 Detection: UV 254 nm
 Temperature: 40 °C
 Sample: Naphthalene
 Inj. volume: 10 μ L

Various organic solvents can be used as RPC mobile phases, but the range of optimal linear velocity varies for different solvents. **Figure 4** shows an H - u curve obtained using a 2.0-mm I.D. column with either methanol or acetonitrile as a mobile phase. Since acetonitrile is less viscous, maximum column efficiency occurs at a higher linear velocity, and high column efficiency is maintained over a greater range of linear velocity. **Figure 5** compares chromatograms obtained using methanol and acetonitrile at flow rates yielding favorable column efficiency (0.20 and 0.50 mL/min, respectively). This figure shows that using acetonitrile, RPC can be performed in about 2/5 of the time taken for methanol without compromising column efficiency. In addition, **Fig. 6** shows the durability of TSKgel ODS-100V 3 μ m column using acetonitrile at a flow rate of 0.50 mL/min. After 1,000 hours, the theoretical plate number had not decreased, and favorable column performance (theoretical plate number and asymmetry factor) was obtained. In this manner, without compromising column efficiency, RPC can be performed in a short time by using a low-viscosity solution as a mobile phase and increasing flow rate.




Figure 4. Effects of linear velocity on height equivalent to a theoretical plate: Effect of organic solvent

Column: TSKgel ODS-100V 3 μ m (2.0 mm I.D. \times 15 cm)
 Eluent: H₂O/CH₃OH = 30/70
 H₂O/CH₃CN = 40/60
 Detection: UV 254 nm
 Temperature: 25 °C
 Sample: Naphthalene
 Inj. volume: 2 μ L




Figure 5. Comparison of chromatograms for standard chemicals

Column: TSKgel ODS-100V 3 μ m (2.0 mmI.D. \times 15 cm)
 Eluent: A) H₂O/CH₃CN = 40/60
 B) H₂O/CH₃OH = 30/70
 Flow rate: A) 0.50 mL/min
 B) 0.20 mL/min
 Detection: UV 254 nm
 Temperature: 25 °C
 Sample: 1. Uracil
 2. Caffeine
 3. Phenol
 4. Methyl benzoate
 5. Benzene
 6. Toluene
 Inj. volume: 2 μ L




Figure 6. Durability under long term flushing with eluent

Column: TSKgel ODS-100V 3 μ m (2.0 mmI.D. \times 15 cm)
 Eluent: H₂O/CH₃CN = 40/60
 Flow rate: 0.50 mL/min
 Detection: UV 254 nm
 Temperature: 25 °C
 Sample: Toluene
 Inj. volume: 2 μ L

2-3 Residual ion exchange activity

In ODS packing materials (C18-bonded silica gel), residual silanol groups affect the retention and peak shape of basic compounds. Similarly to TSKgel ODS-100V 5 μ m, end-capping of residual silanol groups is efficient for TSKgel ODS-100V 3 μ m. Chromatograms for desipramine (basic) and benzene (neutral) before and after end-capping of residual silanol groups in TSKgel ODS-100V 3 μ m were compared (Fig. 7). No marked changes were observed in the retention and peak shape of benzene before and after end-capping (Peak 3). Because desipramine (Peak 2) has electrostatic interactions with residual silanol groups, the chromatogram on TSKgel ODS-100V 3 μ m without end-capping exhibited strong retention and peak tailing (Chromatogram B); in contrast, when TSKgel ODS-100V 3 μ m with end-capping was used, elution and peak shapes were normal (Chromatogram A).




Figure 7. Comparison of chromatograms for basic compounds (desipramine): Effect of end-capping

Column: A) TSKgel ODS-100V 3 μ m (4.6 mm \times 15 cm)
 B) TSKgel ODS-100V 3 μ m (4.6 mm \times 15 cm)
 (Not endcapped)
 Eluent: 5 mmol/L HCOONH₄/ CH₃OH = 20/80
 Flow rate: 1.0 mL/min
 Detection: UV 254 nm
 Temperature: 40 °C
 Sample: 1. Uracyl
 2. Desipramine (52 mg/L)
 3. Benzene
 Inj. volume: 10 μ L

Figure 8 shows the retention of benzene and desipramine using TSKgel ODS-100V 3 μ m with mobile phases of various pH. The retention of benzene, a neutral compound, was mostly stable regardless of mobile phase pH; in contrast, the retention of desipramine, a basic compound, increased with pH due to decreasing amino group dissociation and higher hydrophobicity. Figures 9 and 10 show the retention and asymmetry factors of desipramine using TSKgel ODS-100V 3 μ m and an ODS column with insufficient end-capping at various mobile phase pH values. In general, as shown in Fig. 9, differences in residual silanol group end-capping affect the retention in neutral mobile phases; retention of desipramine was higher for the packing material with insufficient end-capping. Also, as shown in Fig. 10, while use of a packing material with insufficient end-capping in a neutral mobile phase resulted in poor peak shapes (tailing) for desipramine, TSKgel ODS-100V 3 μ m yielded favorable peak shapes with minimal tailing regardless of mobile phase pH.





Figure 8. Relationship between mobile phase pH and retention: Comparison of basic and neutral compounds

Column: TSKgel ODS-100V 3 μ m (4.6 mm \times 15 cm)
 Eluent: 50 mmol/L phosphate buffer (pH 2-7.5)
 /CH₃OH = 30/70
 Flow rate: 1.0 mL/min
 Detection: UV 254 nm
 Temperature: 40 °C
 Sample: Desipramine, Benzene
 Inj. volume: 10 μ L






Figure 9. Relationship between mobile phase pH and retention for a basic compound (desipramine)

Column: TSKgel ODS-100V 3 μ m (4.6 mm I.D.× 15 cm)
 ODS column (4.6 mm I.D.× 15 cm)
 (insufficient end-capping)

Eluent: 50mmol/L phosphate buffer (pH 2-7.5)
 /CH₃OH = 30/70

Flow rate: 1.0 mL/min

Detection: UV 254 nm

Temperature: 40 °C

Sample: Desipramine

Inj. volume: 10 μ L





Figure 10. Relationship between mobile phase pH and peak shape for a basic compound (desipramine)

Column: TSKgel ODS-100V 3 μ m (4.6 mm I.D.× 15 cm)
 ODS column (4.6 mm I.D.× 15 cm)
 (insufficient end-capping)

Eluent: 50 mmol/L phosphate buffer (pH 2-7.5)
 /CH₃OH = 30/70


Flow rate: 1.0 mL/min

Detection: UV 254 nm

Temperature: 40 °C

Sample: Desipramine

Inj. volume: 10 μ L



2-4 Effects of mobile phase on LC/MS(/MS) analysis

In RPC, phosphate buffer is generally used to adjust the pH of the mobile phase. However, in LC/MS(/MS), the use of phosphate buffer, which is nonvolatile, lowers ionization efficiency and contaminates the MS detector; for this reason, it is necessary to use a low-concentration volatile buffer such as formic acid, ammonium formate or ammonium acetate.

Figure 11 shows the effect of mobile phase salt concentration (phosphate and ammonium formate) on the peak shape of desipramine, a basic compound. When using phosphate buffer with high ionic strength as a mobile phase, peak shapes are generally favorable, or the asymmetry factor A_s is close to 1 at buffer concentrations of 5 - 50 mmol/L, independently of buffer concentration, but when ammonium formate, with low ionic strength, was used as a mobile phase, lower buffer concentrations resulted in larger A_s values and greater peak tailing. When using a UV detector in actual testing, favorable peak shapes are obtained by increasing the salt concentration of the mobile phase, but in LC/MS(/MS), the salt concentration of the mobile phase must be low to avoid poor ionization efficiency and the possibility of ion source contamination. In general, the mobile phase concentration is set at ≤ 10 mmol/L.

Figure 12 shows the chromatograms obtained by




Figure 11. Relationship between mobile phase salt concentration and peak shape: Effect of salt type

Column: TSKgel ODS-100V 3 μ m (4.6 mmI.D. \times 15 cm)
 Eluent: 5-50 mmol/L HCOONH₄/CH₃OH = 30/70

5-50 mmol/L Phosphate buffer (pH 7.0)
 / CH₃OH = 30/70


Flow rate: 1.0 mL/min

Detection: UV 254 nm

Temperature: 40 °C

Sample: Desipramine

Inj. volume: 10 μ L



subjecting desipramine (basic compound, Peak 2) to RPC using TSKgel ODS-100V 3 μ m and a competitor's ODS (3 μ m) with 5 mmol/L ammonium formate as a mobile phase. For TSKgel ODS-100V 3 μ m, favorable peak shapes were obtained even when the concentration of ammonium formate was low.

When analyzing a basic compound using low-concentration formic acid or ammonium formate as a mobile phase, the sample concentration has a marked effect on peak shape. **Figure 13** shows the relationships between sample concentration and peak shape (asymmetry factor) obtained by carrying out RPC using TSKgel ODS-100V 3 μ m and phosphate buffer or ammonium formate as a mobile phase. When ammonium formate (low ionic strength) is used, a higher sample concentration (for desipramine) results in greater peak tailing compared to the use of phosphate buffer (high ionic strength). **Figure 14** shows the relationship between sample concentration and peak shape (asymmetry factor) for desipramine using TSKgel ODS-100V 3 μ m and a competitor's ODS (3 μ m), with low-concentration (5 mmol/L) ammonium formate as a mobile phase. When the competitor's ODS (3 μ m) was used, marked peak tailing was seen, starting at the low-concentration area (see **Fig. 12**). As mentioned above, TSKgel ODS-100V 3 μ m exhibits superior properties under the conditions most frequently used in LC/MS(/MS).




Figure 12. Comparison of chromatograms of desipramine obtained using different columns

Column: A) Competitor's ODS (3 μ m)
 (4.6 mmI.D. \times 15 cm)

B) TSKgel ODS-100V 3 μ m (4.6 mmI.D. \times 15 cm)

Eluent: 5 mmol/L HCOONH₄/ CH₃OH = 30/70

Flow rate: 1.0 mL/min

Detection: UV 254 nm

Temperature: 40 °C

Sample: 1. Uracil

2. Desipramine (26 μ g/mL)

3. Benzene

Inj. volume: 10 μ L




Figure 13. Relationship between sample concentration and peak shape: Effect of salt type

Column: TSKgel ODS-100V 3 μ m (4.6 mmI.D. \times 15 cm)
 Eluent: 5 mmol/L HCOONH₄/ CH₃OH = 30/70
 5 mmol/L Phosphate buffer (pH 7.0)
 / CH₃OH =30/70
 Flow rate: 1.0 mL/min
 Detection: UV 254nm
 Temperature: 40 °C
 Sample: Desipramine
 Inj. volume: 10 μ L




Figure 14. Relationship between sample concentration and peak shape: Comparison with competitor's ODS column

Column: A) TSKgel ODS-100V 3 μ m (4.6 mmI.D. \times 15 cm)
 B) Competitor's ODS (3 μ m)
 (4.6 mmI.D. \times 15 cm)
 Eluent: 5 mmol/L HCOONH₄/ CH₃OH = 30/70
 Flow rate: 1.0 mL/min
 Detection: UV 254 nm
 Temperature: 40 °C
 Sample: Desipramine
 Inj. volume: 10 μ L

3. Applications

Figures 15 and 16 are comparison data for TSKgel ODS-100V 3 μ m and TSKgel ODS-100V 5 μ m, while Figs. 17 - 19 show LC/MS chromatograms obtained using TSKgel ODS-100V 3 μ m.

Figures 15 and 16 show chromatograms of acidic and basic compounds measured using TSKgel ODS-100V 3 μ m and TSKgel ODS-100V 5 μ m. Because TSKgel ODS-100V 3 μ m has a higher theoretical plate number, sharper peaks were obtained. Also, the retentions of both samples were mostly comparable regardless of the particle size of the packing material, and no peak tailing was seen for either TSKgel ODS-100V 3 μ m or TSKgel ODS-100V 5 μ m.

Figure 17 shows SIM chromatograms obtained by subjecting aminoglycoside antibiotics to LC/MS using TSKgel ODS-100V 3 μ m. Because aminoglycoside antibiotics are generally highly hydrophilic, sufficient retention cannot be achieved in RPC, so an ion pair reagent is added to the mobile phase. Furthermore, if using an MS detector, the ion pair reagent must be volatile. When we performed LC/MS with addition of heptafluorobutyric acid (HFBA), the peak shapes for all five aminoglycoside antibiotics were favorable.

Figure 18 shows an SIM chromatogram obtained by subjecting microcystin to LC/MS using TSKgel ODS-100V 3 μ m. Microcystin is a hepatotoxin which is synthesized in algal blooms formed due to eutrophication of lakes, and it also acts as a carcinogenic promoter. Therefore, it is necessary to monitor this compound in lake water. The 2003 revision of the drinking water test method requires an assay for microcystin-LR, which is produced by various types of algal bloom. Although the chromatogram was obtained by conducting chromatography at the maximum residue limit (MRL), separation and detection are favorable.

Figure 19 shows SIM chromatograms of sulfonamides obtained by LC/MS using TSKgel ODS-100V 3 μ m. Sulfonamides are synthetic antibiotics which are widely used in veterinary science. The simultaneous analysis methods issued by the Ministry of Health, Labor and Welfare of Japan (Shokuan No. 1129002) "Simultaneous test method (I) for veterinary pharmaceuticals by HPLC (animal husbandry and aquatic products)" lists 16 sulfonamides, and the test method specifies that HPLC be used for quantification and LC/MS(MS) for confirmation. In this experiment, sharp peaks without tailing were obtained for all sulfonamides.




Figure 15. Chromatograms of organic acids

Column: A) TSKgel ODS-100V 3 μ m
(4.6 mm I.D.×15 cm)
B) TSKgel ODS-100V 5 μ m
(4.6 mm I.D.×15 cm)
Eluent: H₂O/CH₃OH = 98/2 + 0.1 % H₃PO₄
Flow rate: 1.0 mL/min
Detection: UV210 nm
Temperature: 40 °C
Sample: 1.Formic acid, 2.Acetic acid
Inj. volume:10 μ L




Figure 16. Chromatograms of basic compounds

Column: A) TSKgel ODS-100V 3 μ m
(4.6 mm I.D. ×15 cm)
B) TSKgel ODS-100V 5 μ m
(4.6 mm I.D. × 15 cm)
Eluent: 50mmol/L phosphate buffer (pH 7.0)
/CH₃OH = 30/70
Flow rate: 1.0 mL/min
Detection: UV254 nm
Temperature: 40 °C
Sample: 1.Uracil, 2.Desipramine, 3.Imipramine
Inj. volume:10 μ L




Figure 17. LC/MS analysis of aminoglycoside antibiotics

Column: TSKgel ODS-100V 3 μ m
(2.0 mmI.D. \times 15 cm)
 Eluent: A: 5mM HFBA
 B: CH₃CN
 Gradient: 0 min (B 10%)
 10 min (B 60%)
 15 min (B 60%)
 Flow rate: 0.2 mL/min
 Detection: MS QTRAP (Applied Biosystems)
 Ion source: ESI
 Polarity: Positive
 Sample: Ribostamycin, Spectinomycin, Gentamicin,
 Bekanamycin, Tobramycin, Kanamycin
 Inj. volume: 5 μ L
 Sample concentration: 0.1 ppm each




Figure 18. LC/MS of microcystin

Column: TSKgel ODS-100V 3 μ m
(2.0 mmI.D. \times 15 cm)
 Eluent: A: 0.1% HCOOH
 B: 0.1% HCOOH in CH₃CN
 Gradient: 0 min (B: 10%)
 10 min (B: 60%)
 15 min (B: 60%)
 Flow rate: 0.2 mL/min
 Detection: MS Q TRAP (Applied Biosystems)
 Ion source: ESI
 Polarity: Positive
 Temperature: 40 $^{\circ}$ C
 Sample: Mycrocystin RR, YR, LR
 Inj. volume: 5 μ L




Figure 19. Simultaneous analysis of sulfonamides by LC/MS

Column: TSKgel ODS-100V 3 μ m (2.0 mmI.D. \times 15 cm)

Eluent: A: 0.1 % HCOOH

B: 0.1 % HCOOH in CH₃OH

Gradient: 0min 0 % B

20min 70 % B

22min 70 % B

23min 0 % B

Flow rate: 0.2 mL/min

Detection: MS Q Trap (Applied Biosystems)

Ion source: ESI

Polarity: positive

Mode: SIM

Temperature: 500°C

Ion spray voltage: 5000V

Temperature: 40°C

Sample: 1. sulfaguanidine 2. sulfacetamide

3. sulfadiazine 4. sulfathiazole

5. sulfapyridine 6. sulfamerazine

7. trimethoprim 8. sulfadimidin

9. sulfamethoxypyridazine

10. sulfamonomethoxine

11. sulfachloropyridazine

12. sulfamethoxazole 13. sulfadoxine

14. sulfadimethoxine 15. sulfabenzamide

16. sulfaquinoxaline 17. sulfanitran

Inj. volume: 2 μ L

(Q TRAP is a registered trademark of Applied Biosystems/MDS SCIEX)

4. Conclusions

As stated above, TSKgel ODS-100V 3 μ m possesses the same separation properties as TSKgel ODS-100V 5 μ m, including high retention of hydrophilic compounds and favorable peak shape for basic compounds. In addition, column efficiency is higher due to the smaller particle size, and high column efficiency can be achieved over a wider range of flow rates, leading to shorter chromatography times. Tosoh has developed various columns with internal diameters ranging from 1.0 to 4.6 mm, including columns for LC/MS(/MS) and microanalysis, allowing selection of the most suitable column size for various applications.