TOSOH

SEPARATION REPORT

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
TSKgel Phenyl-5PW for Use in Hydrophobic Interaction Chromatography of Proteins

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Typical Eluting Condition

Gradient elution of $(\text{NH}_4)_2\text{SO}_4$ concentration from 1 - 2 M to 0 in common buffer solution

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
Change In Phenyl Group Content of TSKgel Phenyl-5PW during the Exposure in Alkaline or Acid Solution for 10 Days at 25°C.

Solution	Phenyl group content (mmol/ml)	
	Before exposure	After exposure
0.5 N NaOH	0.105	0.104
20 % CH_3COOH	0.105	0.106

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
The amount of protein adsorbed on TSKgel Phenyl-5PW in 0.1 M phosphate buffer (pH 7.0) containing various concentrations of ammonium sulfate.

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Dependence of elution volume on ammonium sulfate concentration in initial buffer in high-performance hydrophobic interaction chromatography on TSKgel Phenyl-5PW.

Proteins were separated with 60 min linear gradient of decreasing ammonium sulfate concentration from 1 - 3 M to 0 in 0.1 M phosphate buffer (pH 7.0) at a flow rate of 1 ml/min.

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
Recovery of Enzymatic Activity from TSKgel Phenyl-5PW

Enzyme	Recovery (%)
α -Chymotrypsin (0.4 mg)	92
Lipoxidase (1.0 mg)	89
Ferredoxin NADP reductase (3.0 mg)	100
β -Amylase (1.3 mg)	80
Lysozyme (0.2 mg)	90
Lysozyme (0.05 mg)	90

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
Dependence of peak width on sample loading in high-performance hydrophobic interaction chromatography on TSKgel Phenyl-5PW.

Proteins were separated with 60 min linear gradient of decreasing ammonium sulfate concentration from 1.5 M to 0 in 0.1 M phosphate buffer (pH 7.0) at a flow rate of 0.5 ml/min.

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
Separation of a protein mixture by high-performance hydrophobic interaction chromatography on TSKgel Phenyl-5PW.

Column size: 75 x 7.5 mm I.D.
 Sample: 1. Myoglobin (0.1 or 0.8 mg), 2. Ribonuclease (0.2 mg),
 3. Ovalbumin (0.2 mg), 4. α -Chymotrypsin (0.1 or 0.8 mg)
 Elution: 60 min linear gradient of decreasing ammonium sulfate concentration from 1.5 M to 0 in 0.1 M phosphate buffer (pH 7.0)
 Flow rate: 0.5 ml/min Temperature: 25°C Detection: UV absorbance at 280 nm

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
Separation of a protein mixture by high-performance hydrophobic interaction chromatography on TSKgel Phenyl-5PW

Column size: 75 x 7.5 mm I.D.
 Sample: 1. Cytochrome C, 2. Myoglobin, 3. Ribonuclease, 4. Lysozyme,
 5. α -Chymotrypsinogen, 6. α -Chymotrypsin
 Elution: 60 min linear gradient of $(\text{NH}_4)_2\text{SO}_4$ concentration
 from 1.7 M to 0 in 0.1 M phosphate buffer(pH 7.0)
 Flow rate: 1 ml/min Temperature: 25°C
 Detection: UV absorbance at 280 nm

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
Separation of commercial lipoxidase by high-performance hydrophobic interaction chromatography on TSKgel Phenyl-5PW.

Column size: 75 x 7.5 mm I.D. Sample loading: 1 mg in 0.1 ml
 Elution: 60 min linear gradient of decreasing ammonium sulfate concentration
 from 1.5 M to 0 in 0.1 M phosphate buffer(pH 7.0)
 Flow rate: 0.5 ml/min Temperature: 25°C Detection: UV absorbance at 280 nm
 Recovery of enzymatic activity was 89 %.

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
Separation of commercial α -amylase(barley) by high-performance hydrophobic interaction chromatography on TSKgel Phenyl-5PW.

Column size: 75 x 7.5 mm I.D. Sample loading: 1.3 mg in 0.1 ml
 Elution: 60 min linear gradient from 0.1 M phosphate buffer containing 1.5 M
 ammonium sulfate(pH 7.0) to 0.1 M phosphate buffer(pH 7.0) containing
 5 % 2-propanol
 Flow rate: 1 ml/min Temperature: 25°C Detection: UV absorbance at 230 nm
 Recovery of enzymatic activity was 80 %.

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
Separation of crude ferredoxin NADP reductase by high-performance hydrophobic interaction chromatography on TSKgel Phenyl-5PW.

Column size: 75 x 7.5 mm I.D. Sample loading: 3 mg in 0.1 ml
 Elution: 60 min linear gradient of decreasing ammonium sulfate concentration from
 1.5 M to 0 in 0.1 M phosphate buffer(pH 7.0)
 Flow rate: 0.5 ml/min Temperature: 25°C Detection: UV absorbance at 280 nm
 Recovery of enzymatic activity was 100 %.

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Separation of purified ferredoxin NADP reductase by high-performance hydrophobic interaction chromatography on TSKgel Phenyl-5PW.

Column size: 75 x 7.5 mm I.D. Sample loading: 0.2 mg in 0.1 ml
 Elution: 60 min linear gradient of decreasing ammonium sulfate concentration
 from 1.5 M to 0 in 0.1 M phosphate buffer(pH 7.0)
 Flow rate: 0.5 ml/min Temperature: 25°C Detection: UV absorbance at :

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