SEPARATION REPORT

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- 1 High-performance metal chelate affinity chromatography of proteins. Y.Kato, K.Nakamura and T.Hashimoto, *J.Chromatogr.*, 354(1986)511
- 2 High-performance chromatographic method for the purification of tissue-type plasminogen activator. O.Matsuo, Y.Tanbara, K.Okada, H.Fukao, H.Bando, T.Sakai, *J.Chromatogr.*, 369(1986)391
- 3 Ion-exchange and hydrophobic interaction high-performance liquid chromatography of proteins. D.Josic, W.Hofman and W.Peutter, *J.Chromatogr.*, 371(1986)43
- 4 High-performance liquid chromatography of membrane proteins. Y.Kato, T.Kitamura, K.Nakamura, A.Mitsui, Y.Yamasaki and T.Hashimoto *J.Chromatogr.*, 391(1987)395
- 5 Purification strategies for sendai virus membrane proteins. G.W.Welling, K.Slopsema and S.W.Wester, *J.Chromatogr.*, 397(1978)165
- 6 Interaction of protein with immobilized Cu^{2+} . Quantitation of adsorption capacity, adsorption isotherm and equilibrium constants by frontal analysis. M.Belew, T.-T.Yip, L.Andersson and J.Porath, *J.Chromatogr.*, 403(1987)197
- 7 High-performance analytical applications of immobilized metal ion affinity chromatography. M. Belew, T.T.Yip, L.Andersson and R.Ehrnstrom, *Anal. Biochem.*, 164(1987)457
- 8 High-performance immobilized metal ion affinity chromatography of peptides: Analytical separation of biologically active synthetic peptides. Y. Nalagawa, T.-T.Yip, M.Belew and J.Porath, *Anal. Biochem.*, 168(1988)75
- 9 Characterization of plasminogen activator produced by an established cell line from human ovary, C.Matsuo, H.Fukao and O.Matsuo, *J. Cell. Physiol.*, 134(1988)253
- 10 High-performance immobilized-metal-ion affinity chromatography of peptides and proteins. J. Porath, *J.Chromatogr.*, 443(1988)3
- 11 Evaluation of the interaction of peptides with Cu (II), Ni(II), and Zn(II) by high-performance immobilized metal ion affinity chromatography. T.-T.Yip, Y.Nakagawa and J.Porath, *Anal. Biochem.*, 183(1989)159
- 12 Immobilized-metal affinity and hydroxyapatite chromatography of genetically engineered subtilisin. R.M.Chicz and F.E.Regnier, *Anal. Chem.*, 61(1989)1742