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Ref.; WHO, Draft landscape of COVID-19 candidate vaccines, <https://www.who.int/publications/m/item/draft-landscape-of-covid-19-candidate-vaccines>

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1. J.-H. Tian et al., SARS-CoV-2 spike glycoprotein vaccine candidate NVX-CoV2373 immunogenicity in baboons and protection in mice, Nature Communications (2021)12:372, <https://doi.org/10.1038/s41467-020-20653-8>
2. Q.-D. Su et al., Recombinant SARS-CoV-2 RBD with a built in T helper epitope induces strong neutralization antibody response, Vaccine, 39, (2021) 1241-1247, <https://doi.org/10.1016/j.vaccine.2021.01.044>
3. L. Dai et al., A Universal Design of Betacoronavirus Vaccines against COVID-19, MERS, and SARS, Cell. 2020 Aug 6; 182(3): 722–733.e11., doi: 10.1016/j.cell.2020.06.035, <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC7321023/>
4. W.-H. Chen et al, Optimization of the production process and the characterization of yeast-expressed SARS-CoV recombinant receptor-binding domain (RBD219-N1), a SARS vaccine candidate, J. Pharmaceutical Sciences, 106 (2017) 1961-1970, <https://www.sciencedirect.com/science/article/abs/pii/S0022354917302733?via%3Dihub>
5. W.-H. Chen et al., Yeast-expressed recombinant protein of the receptor-binding domain in SARS-CoV spike protein with deglycosylated forms as a SARS vaccine candidate, Human Vaccine & Immunotherapeutics10:3 648-658 (2014), <https://dx.doi.org/10.4161/hv.27464>
6. European Medicines Agency Assessment report, Covid-19 vaccine, AstraZeneca, EMA/94907/2021, https://www.ema.europa.eu/en/documents/assessment-report/vaxzevria-previously-covid-19-vaccine-astrazeneca-epar-public-assessment-report_en.pdf
7. European Medicines Agency Assessment report, Covid-19 vaccine, Covid-19 Vaccine Janssen, EMA/158424/2021, https://www.ema.europa.eu/en/documents/assessment-report/covid-19-vaccine-janssen-epar-public-assessment-report_en.pdf
8. China Patent CN111218459B (2020); Recombinant novel coronavirus vaccine taking human replication defective adenovirus as carrier, Claim 14. The method of claim 9, wherein the purification method of step (6) is Source 30Q chromatography. <https://www.incopat.com/detail/init2?formerQuery=3eQEo0gaDTqd8BXMTGFaiW4kAd0KKkq&local=en>
9. Sputnik V Press release, April 28, 2021, It uses a 4-stage purification technology that includes two stages of chromatography and two stages of tangential flow filtration. <https://sputnikvaccine.com/newsroom/pressreleases/sputnik-v-statement-on-brazilian-health-regulator-anvisa-s-decision-to-postpone-authorization/>
10. Pharmar's Almanac, April 23, 2020 PAO-M04-20-CL-002, Preliminary scalable downstream processes have been designed, including concentration and buffer exchange of clarified extract capture and polishing chromatography. <https://www.pharmasalmanac.com/articles/applying-rapid-countermeasure-preparedness-to-development-of-a-novel-covid-19-vaccine>
11. Q. Gao et al., Development of an inactivated vaccine for SARS-CoV-2, Science, 03 Jul. 2020: Vol. 369, Issue 6499, pp. 77-81, DOI: 10.1126/science.abc1932. <https://science.sciencemag.org/content/369/6499/77>
12. Wang et al., Cell 182 (2020) 713-721, <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC7275151/pdf/main.pdf>
13. Y.-F. Yao et al., Protective Efficacy of Inactivated Vaccine against SARS-CoV-2 Infection in Mice and Non-Human Primates, Virologica Sinica, (2021), <https://doi.org/10.1007/s12250-021-00376-w>
14. B. Ganneru et al., Evaluation of Safety and Immunogenicity of an Adjuvanted, TH-1 Ske wed, Whole Virion 1 Inactivated SARS-CoV -2 Vaccine, BBV152, <https://www.biorxiv.org/content/10.1101/2020.09.09.285445v2.full>
15. E. Qin et al., Vaccine 24 (2006) 1028-1034, Immunogenicity and protective efficacy in monkeys of purified inactivated Vero-cell SARS vaccine, <https://doi.org/10.1016/j.vaccine.2005.06.038>
16. Z. Kis et al., HOW TO MAKE ENOUGH VACCINE FOR THE WORLD IN ONE YEAR, May 26 (2021), <https://mkus3lurbh3lbtq254fzode-wpengine.netdna-ssl.com/wp-content/uploads/mRNA-vaccine-roadmap-May-26-final.pdf>
17. A. Schmidt et al., Digital twin of mRNA-based SARS-CoV-19 vaccine manufacturing towards autonomous operation for improvements in speed, scale, robustness, flexibility and real-time release testing, Processes 2021, 9, 748. <https://doi.org/10.3390/pr9050748>
18. A. B. Vogel et al., BNT162b vaccines protect rhesus macaques from SARS-CoV-2, Nature 592 (2021) 283, <https://doi.org/10.1038/s41586-021-03275-y>
19. K. S. Corbett et al., SARS-Cov-2 mRNA vaccine design enabled by prototype pathogen preparedness, Nature 586 (2020) 567, <https://www.nature.com/articles/s41586-020-2622-0>
20. E. Prompetchara et al., DNA vaccine candidate encoding SARS-CoV-2 spike proteins elicited potent humoral and Th1 cell-mediated immune response in mice, PLOS ONE, March 22, 2021, <https://doi.org/10.1371/journal.pone.0248007>

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