

## Ordering Information

| Part # | Description  | Matrix  | Housing         | ID (mm) | Length (cm) |  |
|--------|--|---------|-----------------|---------|-------------|--|
| 14449  | TSKgel Boronate-5PW Glass, 10 µm, 100 nm                                 | Polymer | Glass           | 5       | 5           |  |
| 13066  | TSKgel Boronate-5PW, 10 µm, 100 nm                                       | Polymer | Stainless Steel | 7.5     | 7.5         |  |
|        |  |         |                 |         |             |  |
| 14451  | TSKgel Glass Guardgel Kit for 5 mm ID TSKgel Boronate-5PW columns, 20 µm | Polymer | Glass           |         |             |  |
| 13125  | TSKgel Guardgel Kit for 7.5 mm ID TSKgel Boronate-5PW column, 20 µm      | Polymer | Stainless Steel |         |             |  |
|        |  |         |                 |         |             |  |
| 20022  | TSKgel BioAssist Chelate, 10 µm, 100 nm                                  | Polymer | PEEK            | 7.8     | 5           |  |
| 14440  | TSKgel Chelate-5PW Glass, 10 µm, 100 nm                                  | Polymer | Glass           | 5       | 5           |  |
| 08645  | TSKgel Chelate-5PW, 10 µm, 100 nm  | Polymer | Stainless Steel | 7.5     | 7.5         |  |
| 08646  | TSKgel Chelate-5PW, 13 µm, 100 nm  | Polymer | Stainless Steel | 21.5    | 15          |  |
|        |  |         |                 |         |             |  |
| 08647  | TSKgel Guardgel Kit for 7.5 mm ID TSKgel Chelate-5PW column, 20 µm       | Polymer | Stainless Steel |         |             |  |
|        |  |         |                 |         |             |  |
| 14443  | TSKgel Heparin-5PW Glass, 10 µm, 100 nm                                  | Polymer | Glass           | 5       | 5           |  |
| 14444  | TSKgel Heparin-5PW Glass, 10 µm, 100 nm                                  | Polymer | Glass           | 8       | 7.5         |  |
| 13064  | TSKgel Heparin-5PW, 10 µm, 100 nm  | Polymer | Stainless Steel | 7.5     | 7.5         |  |
|        |  |         |                 |         |             |  |
| 14455  | TSKgel Tresyl-5PW, 10 µm, 100 nm   | Polymer | Stainless Steel | 6       | 4           |  |
| 14456  | TSKgel Tresyl-5PW, 10 µm, 100 nm   | Polymer | Stainless Steel | 7.5     | 7.5         |  |
|        |  |         |                 |         |             |  |
| 16208  | TSKgel Tresyl-5PW Guardgel, 10 µm, 2 g                                   | Polymer |                 |         |             |  |

